GC Fuji™ IX GP FAST

RADIOAPAC POSTERIOR GLASS IONOMER RESTORATIVE CEMENT IN CAPSULES

GC Fuji IX GP FAST is a fast-setting version of GC Fuji IX GP CAPSULE.

For use only by a dental professional in the indications for use.

COMPOSITION
Powder: Fluoro-alumino-silicate glass, polycrylic acid, pigment liquid. Water: polycrylic acid, carboxylic acid.

INDICATIONS FOR USE
1. Class I and II restorations in deciduous teeth.
2. Non-load bearing Class I and Class II restorations in permanent teeth.
3. Intermediate restorative and base material for heavy stress situation
In Class I and Class C cavities using sandwich laminate technique.
4. Class V and root surface restorations.
5. Core build-up.

CONTRAINDICATIONS
1. Direct pulp capping
2. In rare cases the product may cause sensitivity in some people.
If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE	Fast	Original
Powder/Liquid Ratio (g/g)	0.40/0.11	0.40/0.11
Mixing Time (sec.)	10"	10"
Time of initial extension after mixing	10"	10"
Working Time (min., sec.)	1'15"	2'00"
Net Setting Time (min., sec.)	2'00"	2'20"
Final Finishing/Commencing Time (37° C. after starting mix) (min., sec.)	3'00"	6'00"

ISO 9917-1 Glass polyalkenoate cement, Restoration

Test conditions: Temperature (23 ±1 °C), Relative humidity (50 ±5%)

1. CAVITY PREPARATION
a) Prepare tooth using standard techniques. Extensive mechanical retention is unnecessary. For pulp capping use calcium hydroxide liners.
b) Apply GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) to the bonding surfaces using a cotton pellet or sponge (Fig. 1).
c) Rinse thoroughly with water. Blot away excess water with a cotton pellet or dry by gently blowing with an air syringe. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).

2. CAPSULE ACTIVATION AND MIXING
a) Before activation, shake the capsule or tap its side on a hard surface to mix the powder (Fig. 2).
b) To activate the capsule, push the plunger until it is fully depressed with the main body and hold it down for 2 seconds (Fig. 3).
1) Ensure the plunger is fully pressed to avoid incorrect mixing ratio of powder and liquid.

3. RESTORATIVE TECHNIQUE
a) Immediately set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000 RPM) (Fig. 4).
Note: When using a rotational mixing device, it is recommended not to use the centrifuge option.
b) Make two clicks in the pre-mix capsule then syringe (Fig. 5). The working time is 1 minute 15 seconds from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
c) Within 10 seconds maximum working time, start to extrude the mixture directly into the preparation (Fig. 6).
e) Form the preliminary contour, and cover with a matrix if required.
f) During the first 3 minutes from start of mix extra care should be taken to avoid moisture contamination or drying-out. In this case this cannot be guaranteed, immediately apply or Fuji Varnish (slow dry) or nGC Fuji Coat LC (light cure).

4. FINISHING
a) Final finishing under water spray using standard techniques can begin approximately 3 minutes after start of mixing (Fig. 7).
b) Apply a final coat of GC G-COA T PLUS (light cure) or GC G-COA T PLUS (UV) alternatively, apply GC Fuji VARNISH (slow dry) or GC Fuji COAT LC (light cure) Refer to the respective instructions for use.
c) Instruct the patient not to apply pressure for 1 hour.

STORAGE
Recommended for optimal performance, store in a cool and dark place (4-25°C) (39.2-77.0°F).

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PAP

<http://www.gwgroup.com>
They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory.gce@gc.dental

Undesired effects - Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM, Ltd
Coopers Court Newport Gagnell
Buckinghamshire
MK16 8JS
United Kingdom

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