








DIRECTA ceramir®

Crown & Bridge QuikCap



DIRECTA AB
Filmskiva 8
SE-194 47 Upplands Väsby
Sweden
Tel. +46 8 506 505 75
www.directadental.com

DIRECTA INC
64 Barnards Road, Unit 3
Newtown, CT 06470
United States
Tel. 800 537 8765
www.directausa.com

CE 2862

01-02-2025-05111

INTENDED PURPOSES

of Ceramir® Crown & Bridge QuikCap

Ceramir® Crown & Bridge devices are radiopaque, bioceramic luting elements intended for permanent cementation of prosthetic teeth or implant abutments in children and adults. The device is used for the fixation of dental crowns and bridges.

Ceramir® Crown & Bridge devices are intended for use by dental professionals for cementation of:

- Metal and porcelain fused to metal crowns and bridges.
- Gold crowns and bridges.
- Cast or prefabricated metal posts.

• High-strength ceramic and bridges suitable for conventional cementation (e.g. zirconia, alumina, and lithium disilicate).

The device is not intended for use as a direct pulp-capping agent.

Users are advised to avoid eye and skin contact, the device contains polyacrylic acid.

Intended users: Preparatory tasks shall be performed by professionally trained dental clinician, dentist or nurse. Application of the cement in the oral cavity of the patient shall be performed by a licensed dental professional.

Intended patient population: Children and adults with teeth:

- Deciduous or permanent teeth.
- Natural or artificial implants (e.g. dental post).

Device description: The cement of Ceramir® Crown & Bridge QuikCap is supplied in a capsule containing cement powder and liquid for mixing and direct application. The capsule is intended for single use.

The cement of Ceramir® Crown & Bridge QuikCap is not intended for reuse after application. The minimum volume of material contained in one capsule is 0.17 ml.

WARNINGS, PRECAUTIONS AND CONTRAINDICATIONS

• Use personal protective equipment, eye protection, face protection, hand protection to avoid contact with the skin, mucous membranes and eyes.

• Do not use on patients with known allergy to polyacrylic acid. In rare cases, the device may cause an allergic reaction. Avoid use of the device in case of known allergy to the material.

• Avoid use of the device without delay after treatment.

• Do not use any pre-treatment agent, e.g. silanes, primers or bonding agent when using the device.

Pre-treatment agents may compromise the bonding of the cement.

• Avoid use of the device in case of eye and skin contact, which may be initiated by eyes, rinses immediately with plenty of water and seek medical care.

• Do not use as a direct pulp-capping agent. In case of dose proximity to the pulp, use a scope and irrigate the pulp to remove the cement.

• In case of skin contact, rinse with plenty of water and in case of contact with eyes, rinse immediately with plenty of water and seek medical care.

• Starting immediately following activation (the reaction between the liquid and powder starts directly after activation).

• Immediately place the capsule in a high frequency oscillating (4.000 to 5.000 rpm) or rotating capsule mixer.

• Immediately insert the capsule into the applicator using the groove to hold it in place.

• Hold the capsule firmly against a hard, flat surface with the plunger facing down (fig. 1).

Push down firmly, forcing the plunger into the body of the capsule until it stops (fig. 1).

• Gently close the dose two times to push the cement towards the tip. Then apply the cement by gently applying the plunger the two more times (fig. 2).

• For a more uniform dispensing, gently tap the applicator button (fig. 3b).

7. Working and setting times

Working time (from end of mixing) at 23 °C/73 °F: 2 minutes.

Setting time (from end of mixing) at 23 °C/73 °F: 2 minutes.

Please observe that room temperatures above 23 °C/73 °F will accelerate setting and reduce working time. Working time can be lengthened by using refrigerated components or by mixing on a cold glass slide.

How to use:

1. Tooth or abutment preparation

Remove the temporary restoration and cement. If needed, clean the tooth using pumice or prophetic powder and/or fluoride.

Rinse with water.

Air dry gently with oil free or cotton rolls. Do not desiccate.

Isolate the site with cotton rolls.

2. Pre-treatment preparation

Lithium Disilicate

Etch with Hf acid according to manufacturer's instructions. Do not use any pre-treatment agents, e.g. silanes, primers or bonding agents.

All other materials: Etch with Hf acid according to manufacturer's instructions.

Do not use any pre-treatment agents, e.g. silanes, primers or bonding agents when cementing with Ceramir® Crown & Bridge. Use of any such pre-treatments will destroy the bond between cement and resin.

An etching agent must be used if recommended by the restoration manufacturer.

3. Pulp protection

If protection is required, use accepted standard clinical procedures.

4. Activation

Start mixing immediately following activation (the reaction between the liquid and powder starts directly after activation).

Open the capsule and remove the cap and remove the capsule.

Hold the capsule firmly against a hard, flat surface with the plunger facing down (fig. 1).

Push down firmly, forcing the plunger into the body of the capsule until it stops (fig. 1).

• Gently close the dose two times to push the cement towards the tip. Then apply the cement by gently applying the plunger the two more times (fig. 2).

• For a more uniform dispensing, gently tap the applicator button (fig. 3b).

8. Sealing and finishing

Seal the margins of the preparation with a rubber-dam or similar material, until the cement reaches a rubber-like consistency (approx. 2 min) (fig. 6).

Remove excess material when the cement has reached the rubber-like consistency (fig. 8).

Check the margins of the preparation for any remaining cement.

Resin composite restorative

Store Ceramir® Crown & Bridge QuikCap capsules at temperatures between 4°C/39°F and 20°C/68°F. Do not open the aluminum foil until immediately before use.

Shelf life:

• Before use and after packaging.

Do not use Ceramir® Crown & Bridge QuikCap after the expiry date.

4. Application

Insert the capsule into the applicator using the groove to hold it in place (fig. 4).

Gently close the dose two times to push the cement towards the tip. Then apply the cement by gently applying the plunger the two more times (fig. 5).

• For a more uniform dispensing, gently tap the applicator button (fig. 3b).

9. Working time

Working time (from end of mixing) at 23 °C/73 °F: 2 minutes.

Please observe that room temperatures above 23 °C/73 °F will accelerate setting and reduce the working time available.

10. Sealing and finishing

Seal the margins of the preparation with a rubber-dam or similar material, until the cement reaches a rubber-like consistency (approx. 2 min) (fig. 6).

Remove excess material when the cement has reached the rubber-like consistency (fig. 8).

Check the margins of the preparation for any remaining cement.

Resin composite restorative

Store Ceramir® Crown & Bridge QuikCap capsules at temperatures between 4°C/39°F and 20°C/68°F. Do not open the aluminum foil until immediately before use.

Shelf life:

• Before use and after packaging.

Do not use Ceramir® Crown & Bridge QuikCap after the expiry date.

4. Application

Insert the capsule into the applicator using the groove to hold it in place (fig. 4).

Gently close the dose two times to push the cement towards the tip. Then apply the cement by gently applying the plunger the two more times (fig. 5).

• For a more uniform dispensing, gently tap the applicator button (fig. 3b).

9. Working time

Working time (from end of mixing) at 23 °C/73 °F: 2 minutes.

Please observe that room temperatures above 23 °C/73 °F will accelerate setting and reduce the working time available.

10. Sealing and finishing

Seal the margins of the preparation with a rubber-dam or similar material, until the cement reaches a rubber-like consistency (approx. 2 min) (fig. 6).

Remove excess material when the cement has reached the rubber-like consistency (fig. 8).

Check the margins of the preparation for any remaining cement.

Resin composite restorative

Store Ceramir® Crown & Bridge QuikCap capsules at temperatures between 4°C/39°F and 20°C/68°F. Do not open the aluminum foil until immediately before use.

Shelf life:

• Before use and after packaging.

Do not use Ceramir® Crown & Bridge QuikCap after the expiry date.

4. Application

Insert the capsule into the applicator using the groove to hold it in place (fig. 4).

Gently close the dose two times to push the cement towards the tip. Then apply the cement by gently applying the plunger the two more times (fig. 5).

• For a more uniform dispensing, gently tap the applicator button (fig. 3b).

9. Working time

Working time (from end of mixing) at 23 °C/73 °F: 2 minutes.

Please observe that room temperatures above 23 °C/73 °F will accelerate setting and reduce the working time available.

10. Sealing and finishing

Seal the margins of the preparation with a rubber-dam or similar material, until the cement reaches a rubber-like consistency (approx. 2 min) (fig. 6).

Remove excess material when the cement has reached the rubber-like consistency (fig. 8).

Check the margins of the preparation for any remaining cement.

Resin composite restorative

Store Ceramir® Crown & Bridge QuikCap capsules at temperatures between 4°C/39°F and 20°C/68°F. Do not open the aluminum foil until immediately before use.

Shelf life:

• Before use and after packaging.

Do not use Ceramir® Crown & Bridge QuikCap after the expiry date.

4. Application

Insert the capsule into the applicator using the groove to hold it in place (fig. 4).

Gently close the dose two times to push the cement towards the tip. Then apply the cement by gently applying the plunger the two more times (fig. 5).

• For a more uniform dispensing, gently tap the applicator button (fig. 3b).

9. Working time

Working time (from end of mixing) at 23 °C/73 °F: 2 minutes.

Please observe that room temperatures above 23 °C/73 °F will accelerate setting and reduce the working time available.

10. Sealing and finishing

Seal the margins of the preparation with a rubber-dam or similar material, until the cement reaches a rubber-like consistency (approx. 2 min) (fig. 6).

Remove excess material when the cement has reached the rubber-like consistency (fig. 8).

Check the margins of the preparation for any remaining cement.

Resin composite restorative

Store Ceramir® Crown & Bridge QuikCap capsules at temperatures between 4°C/39°F and 20°C/68°F. Do not open the aluminum foil until immediately before use.

Shelf life:

• Before use and after packaging.

Do not use Ceramir® Crown & Bridge QuikCap after the expiry date.

4. Application

Insert the capsule into the applicator using the groove to hold it in place (fig. 4).

Gently close the dose two times to push the cement towards the tip. Then apply the cement by gently applying the plunger the two more times (fig. 5).

• For a more uniform dispensing, gently tap the applicator button (fig. 3b).

9. Working time

