Instructions for use

MK-dent Propphy Line Air Polishing Device

Version: 2025-03-19-RevA

Language Index

English Name	Native Name	ISO 639-1	Page
English	English	en	3
French	Français	fr	38
Spanish	Español	es	74
German	Deutsch	de	109
Italian	Italiano	it	144
Chinese	中文	zh	179
Swedish	Svenska	sv	214
Hungarian	Magyar	hu	249
Slovak	Slovenčina	sk	284
Czech	Česky	cs	319
Turkish	Türkçe	tr	354
Norwegian	Norsk	no	389
Bulgarian	Български	bg	424
Lithuanian	Lietuvių	lt	459
Portuguese (Portugal)	Português	pt	494
Arabic	العربية	ar	529
Greek	Αγγλικά	el	564

Table of contents

1	About these instructions for use	6
1.1	Explanation of warnings	6
1.2	Compliance	6
2	Intended purpose.....	7
2.1	Intended use.....	7
2.2	Indications	7
2.3	Contraindications	7
2.4	Intended patient groups.....	7
2.5	Intended users.....	7
2.6	User environment.....	7
3	Safety and requirements	8
3.1	General safety regulations	8
3.2	Reporting obligation	9
3.3	Side effects	9
3.4	Risk patient groups.....	9
4	Product description	10
4.1	Description	10
4.2	Combination with other products	10
4.3	Product components	10
4.4	Scope of delivery	11
4.5	Required materials that are not included	11
4.6	Configurations	11
5	Accessories, consumables and spare parts.....	12
5.1	Accessories	12
5.2	Consumables	12
5.3	Spare parts	12
6	Installation.....	13
6.1	Requirements for the environment	13
6.2	Unpacking.....	13
6.3	Assembly	13
6.4	Preparation for use.....	14
6.4.1	Check and setting of the compressed air.....	15
6.4.2	Connection to the coupling system	15
6.5	Startup and operation	16
6.5.1	Operator protection.....	16

6.5.2	Before use	16
6.5.3	Powder filling	17
6.5.4	Checking for correct operation	18
7	Treatment	19
7.1	Patient, operator, and assistant protection	19
7.2	Spraying	19
7.3	Turning the head	21
7.4	Disconnecting after use	21
8	Processing according to ISO 17664-1	22
8.1	Pre-treatment at the point of use	23
8.1.1	Detaching the head	23
8.1.2	Detaching the powder cap	24
8.1.3	Powder disposal	24
8.1.4	Removal of remaining powder within the air polishing device	24
8.1.5	Draining water	24
8.2	Manual processing	25
8.2.1	Manual cleaning of the powder stirring cap	25
8.2.2	Manual cleaning of the head	25
8.2.3	Manual cleaning of the powder chamber	26
8.2.4	Manual drying	26
8.3	Automatic processing	27
8.3.1	Automatic cleaning, disinfection and rinsing	27
8.3.2	Automatic drying	27
8.4	Sterilization in an autoclave	28
8.4.1	Packaging	28
8.4.2	Sterilization in a steam sterilizer according to EN 13060 / EN ISO 17665	28
8.4.3	Sterilization parameters	29
8.4.4	After Sterilization	29
8.5	Storage	29
9	Maintenance	30
9.1	Regular maintenance and checks	30
9.2	Cleaning in case of problems	30
10	Transport and disposal	30
11	Trouble shooting	31
12	Repair	32
12.1	Replacing the head	32
12.2	Replacing the head O-ring	33

12.3	Replacing the powder cap	33
12.4	Replacing the powder seal	33
13	Technical data	34
14	Explanation of symbols	35
15	Guarantee and warranty	36
15.1	Guarantee	36
15.2	Legal warranty	36
15.3	Conditions	36
15.4	Utilization	37
15.5	Exclusion of liability	37

1 About these instructions for use

These instructions for use apply to all MK-dent air polishing devices (see section **4 Product description**).

All products listed here are only intended for dental treatment in the field of dentistry. Any kind of misuse or modification of the product is not permitted and may lead to hazards. All products listed here may only be used by qualified users in accordance with **2.5 Intended users** for the intended purpose and in compliance with the instructions for use.

1.1 Explanation of warnings

These instructions for use are intended to describe the safe and effective operation of the product to the user. Before operating the product, you must read these instructions for use, observe all safety instructions and warnings and strictly adhere to them. Pay particular attention to all information and procedures described in section **3 Safety and requirements**.

These instructions for use are part of the product and must be kept within reach of the product so that they are always accessible.

To avoid personal injury and damage to property, observe the warnings and safety information contained in this document. These are specially marked:

WARNING

A WARNING alerts you to a potentially serious outcome, unwanted result or safety risk. Failure to observe a warning can result in the death or serious injury of the user or patient.

CAUTION

CAUTION alerts you to situations requiring particular care for the safe and effective use of the product. Disregard for a precautionary measure can result in minor or moderately severe personal injuries or damage to the product or other property and may potentially lead to a remote risk of serious injuries and/or environmental pollution.

NOTE

NOTE

This identifies special tips, such as help to the user or for the improvement of a work process.

1.2 Compliance

This medical device fulfils the requirements of Directive 93/42/EEC.

If you have any further questions regarding the applicable national or international standards, please contact the manufacturer:

MK-dent GmbH, Marie-Curie-Str. 2, 22941 Bargteheide,
Germany
+49 4532 40049-0
info@mk-dent.com

2 Intended purpose

2.1 Intended use

The Prophy Line Air Polishing Device is designed to generate a powder jet with glycine powders that are used during supragingival prophylaxis procedures to remove dental biofilms, plaque, stain and discoloration from tooth surface.

2.2 Indications

The Prophy Line Air Polishing Device is an air polishing device designed for supragingival air polishing as part of the professional cleaning of teeth during prophylaxis treatment by dental professionals to remove dental biofilms, plaque, stain and discoloration from dental enamel.

2.3 Contraindications

Never use the Prophy Line Air Polishing Device on patients with:

- a low salt diet (due to hypernatremia, edema, pregnancy toxemia, etc.),
- a serious ulcer in digestive organs,
- renal dysfunction,
- cardiac dysfunction,
- asthma or chronic pulmonary disease,
- an electrolyte abnormality due to hypochloric alkalosis,
- an intraoral injury or disease,
- hyperemia, inflammation or bleeding in the intraoral cavity,
- oral mucosa inflammation or bleeding,
- allergies
- contact lenses

2.4 Intended patient groups

Patients in need of supragingival air polishing during general prophylactic treatment (professional dental cleaning). However, the decision to use supragingival air polishing is the responsibility of the dentist or dental hygienist.

2.5 Intended users

The Prophy Line Air Polishing Device is an air polishing handpiece for use by dentists or dental hygienists with the appropriate training. Secondary users are the practice staff who are responsible for transport, processing and cleaning.

2.6 User environment

Application takes place in dental practices or dental facilities. Users wear protective gloves. The peripherals used fulfil the requirements of ISO 7494.

3 Safety and requirements

3.1 General safety regulations

WARNING

- This product is supplied non-sterile. Prepare it before first use in accordance with section **8 Processing according to ISO 17664-1**.
- For proper reprocessing, process the device after each patient in accordance with section **8 Processing according to ISO 17664-1**. Safe and effective reprocessing after multiple use without reprocessing the device leads to mixed infections and prevents proper processing due to an atypical accumulation of contamination, occlusions and cleaning agent residues – especially in cavities and supply lines.
- This product must be used only by qualified personnel and only for dental treatment.
- Never disassemble or modify this product.
- This instrument must not be repaired by the user in the event of a malfunction. Contact your dealer or the MK-dent repair service.
- Please do not use the product if the performance of the product deteriorates or in case of malfunctions.
- Do not use this product if it has not been serviced properly.
- Observe the precautionary measures for use in order to ensure correct and safe functioning.
- Use glycine powder for dental prophylaxis only. Other powders (e.g., alumina powder) will cause a decrease in cleaning efficiency, patient injury, powder clogging or damage to the air polishing device.
- Use the powder only once

CAUTION

- The product may only be used in the field of dentistry with a dental treatment unit that has a CE marking or is otherwise certified/approved according to the valid national legal regulations and that is provided with an appropriate water and compressed air supply in compliance with the applicable standards.
- Test for proper functioning before each use and inspect the instrument visually for external damage and contamination.
- Patients with contact lenses must remove them before treatment.
- Depending on the patient's physical condition, allergy symptoms or stomatitis may occur in rare cases. In so, stop using this product immediately.

NOTE

NOTE

Air polishing should be avoided on or near the following:

- Porcelain, polymer, gold alloy and glass ionomer restorations to prevent surface roughness, staining, pitting and the loss of marginal integrity.
- Dental implants, except with glycine or erythritol powder, which can be used safely on implant surfaces.

3.2 Reporting obligation

In the event of serious deterioration in the patient's health during the use of this product, please report this to your specialized dealer or to us as the manufacturer – and to the national competent authority.

Our contact details are included in these instructions for use. The contact details of the notified body will also be provided to you upon request.

3.3 Side effects

Depending on the patient's physical condition, allergy symptoms or stomatitis may seldom occur, depending on the chosen powder. In such a case, stop using this product immediately.

The following side effects may occur if the Air Polishing Device is not used as intended:

- Patient discomfort
- Abrasion of cementum and dentine
- Epithelial abrasion of the gums
- Subcutaneous emphysema

3.4 Risk patient groups

Never use this product on patients with:

- a low salt diet (due to hypernatremia, edema, pregnancy toxemia, etc.),
- a serious ulcer in digestive organs,
- renal dysfunction,
- cardiac dysfunction,
- asthma or chronic pulmonary disease,
- an electrolyte abnormality due to hypochloric alkalosis

Avoid using this product on patients with:

- an intraoral injury or disease,
- hyperemia, inflammation or bleeding in the intraoral cavity,
- oral mucosa inflammation or bleeding,
- allergies,
- contact lenses

4 Product description

4.1 Description

The MK-dent air polishing device is a dental powder jet device to clean tooth surfaces.


- The compact form of the air polishing device handpiece allows easy movement between dental units and takes up little space for storage.
- Thanks to the quick-disconnect system, replacing turbines is easy, enabling immediate use. No complicated installation is required, allowing for simple integration in your treatment system.
- The contra-angled head rotates 360 degrees, giving excellent access to all areas.
- The detachable head makes it easy to rectify powder clogging.
- Driven by compressed air and water from the dental unit, the air polishing device neither affects other medical devices nor is it affected by electromagnetic waves.

4.2 Combination with other products

REF MK-dent air polishing device	PR1011K	PR1011S	PR1011W	PR1011N	PR1011B
MK-dent coupling	QC4012K QC4014K QC5015K QC5016K QC5016KW QC5116K QC5116KW QC6016K QC6016KW QC6116K QC6116KW	QC5016SW QC6016SW	QC4014W QC5016W	QC4014NT QC5016NT QC6016NT	N/A
Connection	KaVo	Sirona	W&H	NSK	BienAir

KaVo, Sirona, W&H, NSK and BienAir are registered trademarks. MK-dent does not have any economical connection to the companies mentioned above.

4.3 Product components


4.4 Scope of delivery

The scope of delivery of the following components and accessories is included in the delivery unit:

- PR1011
- Powder seal (SP3021)
- O-ring head (SP3012)
- Cleaning wire (SP3026)
- Instructions for use (IFU)

When you open the delivery package for the first time, make sure that the seal on the packaging is undamaged and that all components are included in the scope of delivery. If anything is missing, contact your dealer immediately.

4.5 Required materials that are not included

NOTE

NOTE


The scope of delivery does not include powder.
The scope of delivery does not include a coupling.

4.6 Configurations

The configuration of the air polishing devices is as described in **4.3 Product components** – only the connection options differ. The different connections are described as follows:

REF MK-dent air polishing device	PR1011K	PR1011S	PR1011W	PR1011N	PR1011B
Connection	KaVo	Sirona	W&H	NSK	BienAir

KaVo, Sirona, W&H, NSK and BienAir are registered trademarks. MK-dent does not have any economical connection to the companies mentioned above.

5 Accessories, consumables and spare parts

5.1 Accessories

Not applicable

5.2 Consumables

- 3M ESPE, Clinpro Glycine Prophy Powder, grain size 20-25 µm

5.3 Spare parts

CAUTION

The use of unauthorized accessories or unauthorized modifications to the product can lead to injuries or instrument malfunctions.

→ Only use spare parts and components that have been approved by MK-dent for combination with the product.

- Powder seal (SP3021)
- O-ring head (SP3012)
- Cleaning wire (SP3026)

6 Installation

6.1 Requirements for the environment

To supply the products with clean and dry air, clean the compressor air filter in your dental practice. Do not supply moist compressed air. For drainage and regular maintenance of your compressor, please read the operating instructions for your air compressor. Ensure that the compressor is maintained in accordance with its instructions for use.

6.2 Unpacking

When opening the package for the first time, make sure that the content according to **4.4 Scope of delivery** is enclosed with the product. Make sure that the air polishing device has no defects or damage such as cracks.

6.3 Assembly

WARNING

The powder cap and powder seal must be attached firmly. An unstable attachment will cause the powder cap to come off during use and may cause injuries or powder scattering.

CAUTION

Do not yet fill the powder right away. Water from the sterilization process must be expelled by an idle run beforehand.


The air polishing device must be disconnected from the dental unit's turbine hose when attaching and detaching the head. If it is connected, disconnect it from the turbine hose first.

If water, oil or powder adheres to the head or inside of the chamber unit's coupling, wipe it off with a cotton swab or dry cloth.

The powder cap must be attached or detached when the air polishing device is not connected to the turbine hose. If it is connected, detach it from the turbine hose first. If water, oil or powder adheres inside the chamber unit, powder cap or powder seal, wipe it off with a dry cloth.

Step 1:

Attaching the head


Confirm the O-ring of the head is not damaged.

Align the head to the thread of the chamber unit's coupling, turning the head in the direction of the arrow in the illustration until it stops. Confirm that it is attached securely.

Step 2:

Attaching the powder cap

Confirm that the powder seal (blue) is fully engaged in the groove of the chamber unit shown in the illustration.

Turn the powder cap in the direction of the arrow in the illustration. Confirm that it is attached securely.

6.4 Preparation for use

WARNING

Danger of an infection for patients and users due to unsterile instruments
→ Before initial use and after each treatment, the instrument and accessories must be processed according to these instructions for use; see section **8 Processing according to ISO 17664-1**.

Malfunction and/or danger of infection from moist and contaminated compressed air and contaminated water for cooling

→ The product may only be used in the field of dentistry with a dental treatment unit that has a CE marking or is otherwise certified/approved according to the valid national legal regulations and that is provided with a suitable supply of water and compressed in compliance with the applicable standards.


CAUTION

A damaged instrument or damaged components can injure the patient, the user and third parties.

→ Test for proper functioning before each use and inspect the instrument visually for external damage.

→ Do not use the instrument and components if there are any irregularities.

When a handpiece (air turbine, air scaler, etc.) is removed, water and oil will stay behind on the coupling connection. Remove the water and oil completely. Confirm that no water or oil is left behind in the coupling on your turbine hose. If there is any residue left, remove it with a dry cloth.

Make sure that no water or oil remains in the joint of the air polishing device. Remove with a cotton swab if any remains.

6.4.1 Check and setting of the compressed air


WARNING

Excess air pressure will damage dental enamel and cause the patient pain. It may also cause trouble during use such as injuries due to instrument damage or powder scattering.

Connect a pressure gauge to the coupling on your turbine hose until you hear a locking sound. Make sure the measuring instrument is securely attached. If you do not have a pressure gauge, contact your dealer.

Connect the instrument to the pressure gauge coupling until you hear a locking sound. Make sure that the measuring instrument is securely attached.

Turn off the water supply of your dental unit. Direct the air polishing device toward your cuspidor and depress your foot control to activate the dental unit.

Adjust the air pressure of the unit until the gauge shows 0.25 to 0.4 MPa (about 36 to 58 psi; 2.5 to 4.0 bar; 2.5 to 4.1 kgf/cm²). Detach the pressure gauge and the air polishing device after adjusting the pressure.

6.4.2 Connection to the coupling system

WARNING

Malfunction due to decoupling of the instrument during treatment. An instrument that is not properly engaged can come loose from the hose coupling during treatment.

→ Check that the instrument is securely attached to the hose coupling by gently pulling on it before each treatment.

Connect a coupling (sold separately) compatible with your air polishing device to your turbine hose. Confirm that the connection is secure.

6.5 Startup and operation

NOTE

NOTE

Before the start of every working day, the systems conducting water must be flushed for at least two minutes without any instruments attached. Another flushing process must be performed for 20 to 30 seconds after each patient.

6.5.1 Operator protection

WARNING

During use, the operator and assistants must wear gloves, safety goggles and anti-dust masks.

If the powder gets into your eyes, wash your eyes with plenty of water and consult an ophthalmologist.

6.5.2 Before use

CAUTION

Water residue in the channels.

If water remains in the powder chamber, the powder will become moist and cause blockages.

After processing, run the air polishing device for 1 minute without powder before use. This idle run will expel the remaining water inside.

Make sure no powder remains within the powder chamber.

Direct the nozzle toward your cuspidor. Depress your foot control without any water supplied from the dental unit. Run for 1 minute to remove the water remaining in the channels of the air polishing device.

6.5.3 Powder filling

WARNING

Remove powder remaining in/on the screw threads of the powder chamber, powder cap and sealing surface of the powder seal. Remaining powder may cause premature wear of the threads and insufficient sealing, which may cause powder scattering.

If the chamber unit's screw thread fails, including wear and/or deformation, stop using the air polishing device. Request the repair of the chamber unit together with the powder cap.

If the powder cap has a crack or fissure, replace it with a new one.


CAUTION

Do not fill more powder than the prescribed quantity.

Do not shake the air polishing device after filling the powder. The powder may leak from the nozzle.

Always use dry compressed air and dry powder. The use of air and powder containing moisture or oil will cause powder clogging.

Make sure the powder stirring cap is securely in place. Otherwise, the prophylaxis efficiency will be reduced.

- Remove the air polishing device from the turbine hose before (re)filling the powder. Detach the powder cap by turning it in the direction of the arrow in the illustration above.
- Confirm that the powder stirring cap is securely in place.
- If water, oil or powder remains inside of the powder chamber or powder cap, remove it with a dry cloth.
- Fill at most 15 g of powder into the powder cap.
- Take care not to spill the powder over the thread of the powder cap. If this happens, wipe the powder off with a dry cloth.
- Hold the chamber unit downward and the powder cap upward as shown in the illustration. Turn the powder cap in the direction of the arrow as far as it will go. Confirm that the powder cap is screwed securely in place.

6.5.4 Checking for correct operation

WARNING

When connecting your coupling to this product, confirm that the connection is secure after hearing a locking sound. Make sure that the head and powder cap are firmly screwed in. An unstable attachment will risk injuries due to the body and other parts falling off or powder scattering.

Confirm the absence of any damage, cracks, breakage and major corrosion on the product and that it is safe to use.

Never pull the handpiece out of the coupling while the air polishing device is in operation. This can lead to injuries.

- Make sure that the handpiece is firmly in place. Point the nozzle into your drain basin. A large amount of powder may be ejected immediately after filling the powder. Before using it on a patient, spray into your drain basin for a second or two to check for proper function.
- Press the foot switch of your treatment centre once or twice for about one second and check – when air, water and powder are being sprayed – whether the water supply is sufficient and that no leakage of air, powder or water from the powder jet unit can be detected.
- When using the device for the first time, test the spray area and the cleaning effect with an extracted tooth or a tarnished coin.
- If you have problems, read section **11 Troubleshooting**. Use the water volume regulator on your treatment centre to set the water supply.

7 Treatment

7.1 Patient, operator, and assistant protection

WARNING

Before the first contact with the patient, the product must be properly processed to prevent infection. See section **8 Processing according to ISO 17664-1**.

Before using the device, protect the patient with towels and protective eyewear so that the powder will not enter mucous membranes (eyes, nose, tracheae) other than the intraoral mucosa.

During use, the operator and assistants must wear gloves, protective eyewear and anti-dust masks to prevent infection.

If the patient wears contact lenses, the patient must remove them before treatment. If the powder gets into your eyes, wash them with plenty of water and consult an eye doctor.

Protect the patient's face (other than oral cavity) with towels and protective eyewear. When spraying into an area around soft tissues such as gingiva, protect them with cotton in advance.

7.2 Spraying

WARNING

Never spray into subgingival areas, root canals or extraction sockets. This can have undesirable consequences such as emphysema.

Never spray directly into the intraoral soft tissues (gingiva, mucosa), as this can trigger symptoms such as surface damage. Protect the soft tissues with cotton in advance.

Do not put this product close to heat, which would damage the body, risking malfunctioning and possible injuries or powder scattering during use.

CAUTION

Always use the device with sufficient water. Otherwise, the powder will scatter, causing the cleaning efficiency to be reduced.

Keep your vacuum suction system on hand to remove scattered powder. Direct the nozzle toward the tooth surface and depress the foot control to start spraying.


The powder may continue spraying for a few seconds after your foot control is released due to the residual pressure in the powder chamber. Collect the powder with your vacuum suction system to avoid damaging soft tissues. Take care to prevent the patient from swallowing powder.

Avoid prolonged continuous use. Have the patient rinse out his/her mouth as much as necessary.

**CAUTION**

Do not spray directly toward restorations such as crowns, as they may become scratched.

The tooth cuticle is removed during the cleaning procedure, thereby reducing the ability of the tooth to prevent discoloration. The tooth cuticle will be restored by protein contained in saliva. To avoid discoloration, instruct the patient not to smoke or have foods or beverages with a strong coloring effect such as tea, coffee or red wine for two or three hours after treatment.

Keep a distance of 2 to 5 mm from the tooth surface.

Suction the spray with your vacuum suction system and avoid soft tissues.


Spray from the gingival margin to the occlusal surface.

When cleaning a large area, spray while moving the nozzle in minute increments.


The powder may continue spraying for a few seconds after the foot control is released due to the residual pressure within the powder chamber.

Collect the powder with your vacuum suction system to avoid damaging soft tissues.

Have the patient rinse out his/her mouth as needed or after treatment.

Never spray into subgingival areas, root canals or extraction sockets. This can have undesirable consequences such as emphysema.


Never spray directly into the intraoral soft tissues (gingiva, mucosa), as this can trigger symptoms such as surface damage.

7.3 Turning the head

WARNING

Be sure to rotate the head by holding and turning the grip, not the head itself. Otherwise, the head will loosen, risking injuries or powder scattering during use.

Rotate the head by turning the grip.
The head rotates 360 degrees.

7.4 Disconnecting after use


WARNING

The product must not be dropped.

CAUTION

Remove the remaining water, oil or powder in and on the head, the coupling of the chamber unit, the powder chamber and the powder cap. Otherwise, the powder will get damp and cause clogging.

Detach the handpiece according to the coupling's instructions for use.

When reconnecting it, clean the coupling, turbine hose and inside of the chamber unit's joint in advance according to **6.4.3 Connection to the coupling system**.

8 Processing according to ISO 17664-1

NOTE

NOTE

The possible number of applications does not depend on the processing cycles but rather on the use of the device and the resulting wear and tear. The product's service life is therefore determined by signs of wear and tear and damage occurring during use.

Dispose of damaged products in accordance with section **10 Transport and disposal**. Do not use any products with obvious signs of wear or damage.

Only use cleaning agents / disinfectants that are authorized for the processing of medical devices. Follow the manufacturer's instructions for the cleaning agents / disinfectants.

If unsuitable cleaning agents or disinfectants are used, this can have negative effects on the instruments:


- Damage or corrosion
- Discoloration of the product
- Corrosion of metal parts
- Reduced service life

Follow the instructions of the cleaning agent / disinfectant manufacturer. Do not use any abrasive cleaning agents / disinfectants. The manufacturer recommends using only EN ISO 15883-1 and -2-compliant washer-disinfectors for automated cleaning/disinfection.

The head of the instruments can be sterilized. Do not use sterilization temperatures above 135°C (275 °F).

Ultrasonic cleaning is not recommended for processing the instruments, as this could reduce the product's service life.

Assembly and disassembly are described in sections **6.3 Assembly, 8.1.1 Detaching the head** and **8.1.2 Detaching the powder cap**.

8.1 Pre-treatment at the point of use

WARNING

Risk of infection from non-sterile instruments. Always wear protective gloves.

CAUTION

Do not use ultrasonic cleaners to clean this product. This will lead to malfunctions or damage.

Using disinfecting baths and disinfectants containing chlorine can cause defects and malfunctions of the instrument.

The cleaning result is diminished by surface drying of residues.

→ Rinse the instrument immediately after finishing work on the patient to remove residue.

Do not submerge the instrument in solutions or similar.

Remove the powder that remains inside and on the surface of the air polishing device between treatment and sterilization.

NOTE

NOTE

Clean the coupling on your turbine hose and inside of the instrument connector according to **6.4 Preparation for use** before reconnecting the instrument.

- Remove saliva, blood or powdery deposits immediately.
- Do not place the instrument in any kind of liquid.
- Flush the inner tubes of the product according to the RKI specifications for 20 seconds while it is still connected to the treatment unit.
- The instrument should be pre-cleaned directly after the treatment for the sake of personal safety. The aim is to prevent organic material and chemical residue from drying in the cavity or on the outer parts of the instruments and to avoid contaminating the environment.
- Remove adhering contaminants, body fluids and tissue using a disposable cloth / paper tissue.
- Flush the product with water after use (temperature below 35 °C (95 °F)).
- Do not use any abrasive cleaning agents or water with temperatures above 35 °C (95 °F), as this can lead to increasing adhesion of the contaminants, making further processing steps more difficult.
- Clean the instrument within one hour after each treatment to prevent surface drying.
- Ensure that the instrument is stored safely and transported in dry condition in a closed container to the place of processing.

8.1.1 Detaching the head

The head must be attached or detached when the air polishing device is not connected to a turbine hose. If the air polishing device is connected to a turbine hose, disconnect it from the turbine hose. Turn the head in the opposite direction of the arrow as shown in **6.3 Assembly**.

Before reattaching it, clean the head and inside of the chamber unit's coupling in advance according to **6.3 Assembly**.

8.1.2 Detaching the powder cap

The powder cap must be attached or detached when the air polishing device is not connected to a turbine hose. If the air polishing device is connected to a turbine hose, disconnect it from the turbine hose.

Hold the air polishing device, keeping the powder cap downward. Turn the powder cap to detach it. Detach the powder cap above your cuspidor, as the powder in the powder chamber may scatter when the powder cap is released. Discard any leftover powder in the powder cap.

Before refilling the powder and reattaching the powder cap, clean the inside of the powder chamber and powder cap according to **6.3 Assembly**.

8.1.3 Powder disposal


CAUTION

The leftover powder must be discarded after treatment. If the powder gets damp, it will cause clogging.

Detach the powder cap according to **8.1.2 Detaching the powder cap**. Discard any leftover powder.

8.1.4 Removal of remaining powder within the air polishing device

Attach the emptied powder cap to the chamber unit.

Connect the air polishing device to the coupling on your turbine hose.


Direct the nozzle toward your cuspidor.

Depress your foot control without any water supplied from the dental unit.

Run for 1 minute to remove the powder remaining within the channels of the air polishing device.

Disconnect the air polishing device from your coupling.

8.1.5 Draining water

Attach the powder cap and head to the emptied chamber unit.

Direct the nozzle toward your cuspidor.

Blow air into the joint of the chamber unit to drain the water channels.

8.2 Manual processing

CAUTION

Insufficient cleaning of the surface and of water and air supply lines

Thorough manual pre-cleaning must be carried out to ensure successful cleaning.


The use of manual pre-cleaning requires training of the employees who carry out these processing instructions. This is necessary to guarantee successful cleaning.

Malfunctions can occur if cleaning is done in the ultrasonic machine.

→ Clean the instrument exclusively as described below

Cleaning must be done when the air polishing device is disconnected from the coupling on your turbine hose. After disassembling the head and powder cap from the chamber unit wash each part (in particular, where powder adheres) with running water.

8.2.1 Manual cleaning of the powder stirring cap


Remove the powder stirring cap from the powder pipe.

(The cleaning wire does not pass through the center hole when the powder stirring cap is connected.)

Use a cotton swab or dry cloth to remove the powder adhered to the center hole of the powder stirring cap.

Clean the center hole with the attached cleaning wire to remove any powder blocking the center hole.

8.2.2 Manual cleaning of the head

Remove the head from the chamber unit.

Insert the cleaning wire into the spraying orifice of the nozzle and center hole of the connection.

Use the cleaning wire to remove any powder remaining within the channels.

Blow air into the spraying orifice of the nozzle to remove powder remaining inside.

If any water, oil or powder remains at the surface of the head, remove it with a dry cloth.


8.2.3 Manual cleaning of the powder chamber

Detach the head and powder cap from the chamber unit.

Loosen the powder stirring cap from the powder pipe by turning it in the direction of the arrow in the illustration and detach it.

Insert the cleaning wire into the powder pipe and remove the powder remaining in the pipe.

Blow air into the powder pipe of the nozzle to remove any powder remaining inside.

If any water, oil or powder remains inside the chamber unit or in the powder cap, remove it with a cotton swab or dry cloth.

Attach the powder stirring cap to the powder pipe by turning it in the opposite direction of the arrow in the illustration.

CAUTION

Make sure that the powder stirring cap is securely in place after cleaning. Improper attachment will reduce the cleaning efficiency.

8.2.4 Manual drying

Blow on the instrument with compressed air from the inside and the outside until no visible drops are left.

8.3 Automatic processing

8.3.1 Automatic cleaning, disinfection and rinsing

NOTE

NOTE

MK-dent recommends thermal disinfection devices that satisfy the requirements of EN ISO 15883-1 and are approved by the manufacturer for dental handpieces and turbines. The cleaning should be performed at min. 55 °C (131 °F) for at least 6 minutes and disinfection at min. 90 °C (194 °F) for at least 5 minutes (for A0 value > 3000). A mildly alkaline cleaning agent with pH 9 to 11, such as Dr. Weigert Neodisher® MediClean forte, is recommended.

CAUTION

For automatic cleaning of the interior, use the adapters specified by the device manufacturer.

Only process instruments that are dry and free of residue using an automatic cleaning and disinfection process.

The air polishing device can be automatically processed when assembled. For automatic processing, the adapters specified by the manufacturer of the thermal disinfectant and compatible with your air polishing device must be used. For information about the positioning, refer to your device's instructions for use.

The automatic cleaning and disinfection process was validated using the device WD BHT INNOVA® M3 in program no. 03 (Thermal disinfection) with Dr. Weigert Neodisher® MediClean forte 0.5% as the cleaning agent. The Instruments have been connected to the machine using a silicon hose with a Luer-Lock adapter.

Parameters for automatic cleaning and disinfection:

- Pre-rinse for 1 minute with 0.5% Neodisher® MediClean forte detergent at 30 °C (86 °F)
- Clean for 6 minutes with 0.5% Neodisher® MediClean forte detergent at 55 °C (131 °F)
- Rinse for 1 minute with demineralized water
- Thermal disinfection for 5 minutes at 90 °C (194 °F) using low-germ (< 10 CFU/ml) and low-endotoxin (< 0.25 endotoxin units/ml), demineralized water (for an A0 value > 3,000).

8.3.2 Automatic drying

The drying process is normally part of the cleaning process of your thermal disinfectant. Please follow the instructions for your thermal disinfectant. Allow the instrument to cool down before removing it to ensure safe and easier removal from the adapter.

If there is still visible contamination after treatment in the thermal disinfectant, the process must be repeated.

To avoid any kind of damage to the instrument, make sure after each cycle that the instrument is dry inside and outside and free of residue. If moisture remains on or in the instrument after the cleaning process, dry it with a lint-free cloth. Also check the cavities of the instruments for dryness and dry them with medical compressed air if necessary.

8.4 Sterilization in an autoclave

8.4.1 Packaging

NOTE

NOTE

Use a validated packaging procedure that is suitable for sterilization. The sterilization bag must be big enough for the instrument so that the packaging is not stretched. The sterilization packaging must comply with EN ISO 11607-1 and be suitable for the sterilization process. Pack each instrument individually in a simple sterilization bag. Seal the sterilization bag carefully and check the seal for leaks.

Remove the head from the chamber unit. Confirm that no powder remains in the head. Put the head into a single sterilization bag and seal it.

8.4.2 Sterilization in a steam sterilizer according to EN 13060 / EN ISO 17665

WARNING

Only the head of the air polishing device can be sterilized in an autoclave. Sterilization must be carried out before the first use and after each patient to prevent direct and cross-infection.

To avoid infection, wear suitable protective clothing such as a mask, gloves and safety goggles.

Contact corrosion due to moisture can cause damage to the product.

→ Remove the instrument from the steam sterilizer immediately after the sterilization cycle.

This medical device has a maximum temperature resistance of 135 °C (275 °F).

CAUTION

To avoid problems and damage to the product, follow the instructions when sterilizing.

Do not sterilize with any powder remaining.

Do not sterilize the head with instruments that have any drug or chemical agent residue.

Do not immerse the air polishing device in disinfectant liquids.

Sterilizers with dry heat generate high temperatures that exceed the maximum sterilization temperature of this product. Do not use these in order to avoid mechanical problems or damage.

Depending on the conditions of use of your sterilizer or your sterilization methods, this product may undergo discoloration.

NOTE

NOTE

Pay close attention to the manufacturer's instructions for use.

Put the sealed head in the sterilization bag on a sterilization tray to avoid direct contact with the inner chamber surface of the sterilizer and then place it in the autoclave. MK-dent recommends the use of a sterilizer according to EN 13060, EN ISO 17665 or EN 285.

8.4.3 Sterilization parameters

Only use the following sterilization parameters to achieve an SAL of 10^{-6} :

Autoclave with a fractionated pre-vacuum: at least 5 minutes at a minimum of 134 °C (273 °F)

The sterilization was validated with the device W&H LISA 517 in a fractionated pre-vacuum process at 134 °C (273 °F) for 5 minutes. According to the sterilizer manufacturer's instructions for use, the drying time for each sterilization process is calculated by the autoclave's the software.

The instructions listed above for processing a medical device for reuse have been validated and found to be suitable by the medical device manufacturer. The processor is responsible for ensuring that the processing actually performed, including the equipment, materials and personnel used in the processing facility will achieve the intended result. This requires verification and/or validation and routine monitoring of the process.

8.4.4 After Sterilization

When the sterilization cycle is finished, take the sterilized head out of the autoclave. When doing so, be sure to use appropriate personal protective equipment, as the instruments will be hot.

8.5 Storage

After sterilization and drying, check that there are no drops of residual water. Store the processed instruments in a dry, dark and cool room, protected from dust and with the lowest possible microbial load.

Temperature	-20°C to +70°C (-4°F to +158°F)
Relative humidity	5% to 95%, non-condensing
Air pressure	700 hPa to 1060 hPa (10 psi to 15 psi)

CAUTION

The chamber unit and the head must be completely dry before storage. If any water remains, discoloration, rust or powder clogging will occur.

Instrument malfunction during operation after heavily cooled storage conditions
→ Allow heavily cooled instruments to warm up to room temperature before starting to use them.

Do not store the instrument in a place where it could be exposed to water.

Do not store the instrument in places where chemicals are stored or where corrosive gases exist.

9 Maintenance

9.1 Regular maintenance and checks

WARNING

Inappropriate maintenance or service can result in malfunctions and/or deterioration.

→ Always conduct proper maintenance and servicing.

CAUTION

For safe use of the product, perform checks and maintenance of the instrument according to these instructions for use.

9.2 Cleaning in case of problems

In case of problems with the air or water supply lines, clean each part of the instrument. The instrument must be disconnected from the hose coupling for cleaning. Instructions can be found in section **8.1 Pre-treatment at the point of use**.

10 Transport and disposal

CAUTION

Users or third parties may be infected by contaminated medical devices. Process the product before transport or disposal or dispose of it properly in compliance with the applicable regulations for the disposal of contaminated medical waste.

Process the instrument in accordance with section **8 Processing according to ISO 17664-1** and return it to us in the original packaging. We will take care of the professional disposal of the device and the original packaging for you.

11 Trouble shooting

NOTE

NOTE

This instrument must not be repaired by the user in the event of a malfunction.
Please contact your dealer or the MK-dent repair service.

If this instrument does not function as described in these instructions for use, check the list below for your problem. If your problem is not resolved or if your problem is not listed below, please contact your dealer.

Symptoms	Possible causes	Countermeasures
1) No powder spray, low cleaning efficiency	• The center hole of the powder stirring cap is clogged with powder • Powder clogging within the channels • The powder stirring cap is not attached properly • No air supply • Low air pressure at the coupling. • The control system of the dental unit does not work	• Clean the powder stirring cap according to 8.2.1 Manual cleaning powder stirring cap • Clean the head and chamber unit according to 8.2.2 Manual cleaning head and 8.2.3 Manual cleaning of the powder chamber • Attach the powder stirring cap properly according to 8.2.3 Manual cleaning of the powder chamber • Set coupling pressure between 0.25 and 0.4 MPa • Check the control system using another handpiece
2) No water	• The water supply system of your dental unit is shut off • Water does not reach turbine hose • Water channels of air polishing device are clogged • O-ring(s) on your coupling is/are damaged	• Activate the water supply system of your dental unit • Check the water supply system of your dental unit • Ask your dealer for repair • Replace the O-ring(s) on your coupling (sold separately or included in delivery)
3) Leak (air, water or powder) from powder chamber coupling	• The head is not firmly attached • The O-ring in the head is worn or damaged	• Attach the head firmly according to 6.3 Assembly • Replace the attached head O-ring according to 12.2 Replacement of the head O-ring
4) Air and/or powder leak from near the powder chamber	• The powder cap is not firmly attached • The powder seal is not properly attached in the powder chamber • The powder seal is worn or damaged • Powder adheres to the sealing surface (of the powder seal and/or powder cap)	• Attach the powder cap firmly according to 6.3 Assembly • Attach the powder seal properly according to 12.4 Replacement of the powder seal • Replace with the attached new powder seal according to 12.4 Replacement of the powder seal

Symptoms	Possible causes	Countermeasures
		Fill the powder into the powder cap properly according to 6.5.3 Powder filling
5) Water leak from around the coupling on the turbine hose	- O-ring(s) on your coupling is/are damaged (for PR1011K/W/S/B) - O-ring(s) on your joint is/are damaged (for PR1011N)	- Replace the O-ring(s) on your coupling (sold separately or included in the delivery package) - Replace the O-ring(s) on the joint (separately sold)
6) Air polishing device can not be connected to your coupling	- The air polishing device is not compatible with your coupling - The light cover at the front of the coupling is loose (if you are using a coupling with a light)	- Use the air polishing device compatible with your coupling according to 4.2 Combination with other products - Turn the light cover at the front of the coupling to fix it securely

12 Repair


CAUTION

Contaminated medical devices can infect users or third parties. Process the product before dispatching it.

Process this instrument in accordance with section **8 Processing according to ISO 17664-1** and send it to your dealer or directly to MK-dent if it is defective. Please contact MK-dent customer service to arrange for shipping or collection for factory repair. Questions about noise development, maintenance, preparation and functionality can often be clarified in a conversation, making it unnecessary to send in the product.

12.1 Replacing the head


The nozzle of the head will wear out with use. A broadened spraying orifice will reduce the cleaning efficiency and make pinpoint spraying more difficult.

When the spraying orifice has become larger as shown in the illustration, restore the work efficiency by replacing the head (sold separately).

The entire head must be replaced, not only the nozzle. Attach the head properly according to **6.3 Assembly**.

12.2 Replacing the head O-ring

If the O-ring at the coupling of the head is damaged or worn or if air, water or powder is leaking from the joint, replace the used O-ring with the head O-ring provided (spare part).

Head o-ring (spare part)

Secure the O-ring properly in place. Do not damage the surface of the O-ring when replacing it.

12.3 Replacing the powder cap

The transparency of the powder cap will decrease with use and may result in low visibility inside the powder cap. If you cannot see powder in the powder chamber, replace the powder cap (sold separately). If the powder cap is cracked, it must be replaced with a new one. Secure the powder cap properly according to **6.3 Assembly** when replacing.

12.4 Replacing the powder seal

WARNING

When replacing the O-rings and powder seal, ensure that they are fitted in the correct position and direction. Incorrect assembly can lead to injuries due to separation of the head or powder cap, powder scattering and leakage of air, water or powder during use.


If the powder seal fixed in the thread of the powder chamber is damaged or worn or if air or powder leaks between the powder cap and powder seal, replace the powder seal in use with the powder seal provided (spare part). Proceed as follows:

If powder or anything else remains in the groove of the powder chamber where the powder seal is attached (see illustration below), remove it with air from a three-way syringe. Clean the groove and secure the powder seal properly so that the powder seal is fully seated in the groove (see illustration below). Do not damage the surface of the powder seal when replacing it.

Leak

Insert the seal completely into the groove of the chamber unit

Enlarged view

Groove

13 Technical data

Article no.	PR1011K	PR1011S	PR1011W	PR1011N	PR1011B
Device length	201 mm (7.9 inch)	196 mm (7.7 inch)	200 mm (7.9 inch)	208 mm (8.2 inch)	196 mm (7.7 inch)
Weight	142 g (5.0 oz.)	137 g (4.8 oz.)	136 g (4.8 oz.)	144 g (5.1 oz.)	133 g (4.7 oz.)
Connection	KaVo	Sirona	W&H	NSK	BienAir
Powder consumption volume	4 g / min (air pressure on coupler 0.4 MPa / 4 bar)				
Coupling air pressure	0.25 - 0.4 MPa (36 - 58 psi, 2.5 - 4.0 bar, 2.5- 4.1 kgf/cm ²)				
Air consumption volume	10 NI / min (air pressure on coupler 0.4 MPa / 4 bar)				
Water flow volume	40 ml / min (water pressure on coolant 0.2 MPa / 2 bar)				

KaVo, Sirona, W&H, NSK and BienAir are registered trademarks. MK-dent does not have any economical connection to the companies mentioned above.

14 Explanation of symbols

	Caution – see section 1 Explanation of the warnings
	Sterilizable in a steam sterilizer (autoclave) up to the specified temperature
	Suitable for automatic processing in thermal disinfection devices
	Consult instructions for use and refer to the source of the electronic instructions for use
	CE marking with identification number of the notified body – indicates that the product is in compliance with European legislation for medical devices
	Medical device (CE conformity assessment procedure)
	Catalog number
	Serial number
	Manufacturer
	Date of manufacture
	Unique product identification - UDI (Unique Device Identifier) – indicates a barcode as containing Unique Device Identification
	HIBC Code (Health Industry Bar Code)

15 Guarantee and warranty

MK-dent assures you that your instrument has undergone a thorough **quality inspection** before leaving the factory. These quality inspections are carried out according to considerably stricter standards than the requirements of the product standard (ISO 20608 Dentistry – Powder jet handpieces and powders). In addition, we offer you comprehensive warranty services to ensure your satisfaction.

15.1 Guarantee

MK-dent grants the end customer a guarantee of perfect function and freedom from defects in material and workmanship from the date of purchase for a period of 24 months.

Exclusions:

- Wear parts, e.g., collets, ball bearings, gaskets, adapters, bulbs and O-rings
- Fall damage
- Willful damage and improper use contrary to the instructions for use
- Transport damage due to improper packaging (Instruments must be packed individually and with sufficient padding)
- No enclosure of the invoice copy of the device purchase (invoice recipient must be a dental practice)
- The warranty is valid from the date of purchase and is non-transferable
- Defects or their consequences due to interventions or modifications to the product by the end user or by a third party not authorized by MK-dent

15.2 Legal warranty

MK-dent grants the statutory warranty on all products to the contractual partner. This covers material and processing defects that were already present at the time of the transfer of risk or that occur within the statutory warranty period. Warranty claims will only be recognized if the product is submitted with proof of purchase in the form of an invoice or a copy of the delivery note (contractual partner, date of purchase, model and serial number must be clearly visible).

Exclusions:

- Wear parts, e.g., collets, ball bearings, gaskets, adapters, bulbs and O-rings
- Fall damage
- Willful damage and improper use contrary to the instructions for use
- Transport damage due to improper packaging (instruments must be packed individually and with sufficient padding)
- The warranty is valid from the date of purchase and is non-transferable
- Defects or their consequences due to interventions or modifications to the product by the end user or by a third party not authorized by MK-dent

15.3 Conditions

- Proof of purchase (invoice or delivery note) is required for guarantee or warranty claims.
- The guarantee or warranty is void if the product is tampered with or modified by unauthorized third parties.
- Repairs or replacements under warranty do not entitle the warranty period to be extended or restarted.
- MK-dent reserves the right to repair or replace defective instruments.

15.4 Utilization

In the event of a warranty claim, please contact MK-dent customer service directly. In the event of a warranty claim, please contact your contractual partner.

Contact MK-dent customer service

Phone: +49 4532 40049-0

Email: service@mk-dent.com

Please have the following information ready:

- Product model and serial number
- Proof of purchase with date
- Detailed description of the problem

Our customer service team will guide you through the warranty process and ensure a quick and satisfactory solution. After the repair, your instrument will be returned to you securely packed and sealed together with a detailed test report.

15.5 Exclusion of liability

Within the scope of the guarantee or warranty, MK-dent is not liable for financial losses, downtimes, loan or rental equipment, loss of profit or similar consequences. The liability of MK-dent is limited to the purchase value of the product.

MK-dent is not liable for any defects and their consequences that may arise due to natural wear and tear, improper handling, improper cleaning or maintenance, non-compliance with operating or handling instructions, calcination or corrosion, contaminated air and water supply or chemical or electrical influences that are unusual or not permitted according to the instructions for use of MK-dent or other manufacturers' instructions for use.

Any liability is excluded if defects or their consequences are caused by interventions or modifications to the product by the end user or by a third party not authorized by MK-dent.

MK-dent reserves the right to change the warranty conditions at any time. Changes have no retroactive validity for products already purchased.

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