SWAN+ Curing Light Instruction Manual

Please read this manual before operating the device

Guilin Refine Medical Instrument Co., LTD.
SE1-SWP-M001 Version: 1.0 20260112

Safety Precautions

1. This product is intended for use in hospitals and dental clinics only. The use of the product must comply with the requirements of the relevant operating specifications and relevant regulations of the medical department. Users must be professionally trained and qualified as dentists or technicians. Do not look directly into the light output. Patient, clinician and assistants are always suggested wearing protectors (safety goggles, face shield...) to avoid crossing infection. Please do not use any voltage other than the specified. Before connecting the built-in ultrasonic scaler without transformer to power supply, please check the output voltage is 24VAC, in case of connecting to wrong power supply which may break the unit.
2. The products are reusable:
Prior to each use, place a new disposable sleeve over the curing light Swan+. The hygienic disposable sleeve is custom fitted to the curing light and keeps the surface of the curing light clean. The disposable sleeve helps prevent cross contamination, helps keep dental composite material from adhering to the surface of the lens and curing light, and prevents discoloration and corrosion from cleaning solutions. The curing lights must be cleaned and sanitized with appropriate cleaning and/or sanitizing agents after each patient. The sleeve is equipped with the curing light.
3. During the operation, the light should be aimed directly at the composite resin to ensure the effect of solidification.
4. Please use the original light hood and follow the right installation ways to avoid the damage of blue light to the eyes.
5. Only the specified voltage must be used to recharge the battery. Only the original pedestal charger, adapter and lithium battery must be used. Other brand pedestal charger, adapter and lithium battery are likely to damage the circuit. Don't use an impure water source, and be sure not to use normal brine instead of pure water source.
6. Prior to initial use, the battery must be fully charged for a minimum of 4 hours. Continuous connection of the charger to the charging base without active charging is prohibited. Operation of the curing light during the charging process is not permitted.
7. Contact between the charging connector and metal objects or other conductive materials is strictly prohibited to prevent damage to the charging circuit or the battery.
8. Recharge the battery in cool and ventilated room. Press the main unit to the buckle of charging base, in order to avoid the poor connection.
9. Disassembling the battery is strictly prohibited, as it may result in short-circuiting or leakage.
10. Do not crush, shake, or subject the battery to impact. The lithium battery must not be short-circuited or stored in contact with metal objects or other conductive materials.
11. The instrument has electromagnetic interference. Do not use it around electronic surgery, and use the instrument with caution in an environment with strong electromagnetic interference.
12. This product does not contain toxic or hazardous substances and it is discarded in accordance with the relevant laws and regulations on waste medical devices. Please deal with the replaced battery according to national regulations.
13. As a professional manufacturer of medical instruments, we are only responsible for the safety on the following conditions:
 - The maintenance, repair and modification are made by the manufacturer or the authorized dealer.
 - The changed components are original of our company and operated correctly according to instruction manual.
14. The device must not be used in MRI environment for the device is easily affected by the electromagnetic emission and would not work or work normally.
15. When you meet circumstances where the main unit expires, the misuses lead to the short circuit of the circuit board or accidentally dropping the device results in the damage of components, the device should no longer be reused.
16. Special precautions/considerations should be taken before use on patients with heart disease, pregnant woman, children and persons with sensitivities to blue light.

⚠ Please report any serious incident that may occur in connection with this product and its accessories to the manufacturer as well as to your national authority.

Symbol instruction

Symbol	Instruction	Symbol	Instruction
	Caution		Refer to instruction manual/booklet
	Date of manufacture		Manufacturer
	Class II equipment		Type B applied part
	Keep dry		Recovery
	Fragile, handle with care		For indoor use only
	Atmospheric pressure limitation: 70kPa-106kPa		Humidity limitation: 10%-93%
	Temperature limit: -20°C - +55°C		Waste electrical and electronic equipment
	Separate collection of batteries		Importer
	Serial number		Unique device identifier

	Medical device		CE marking with identification number of the Notified Body
	Cross-contamination		Protective glasses and gloves
	Non-sterile Condition		Important Information
	Authorized representative in the European Community		Power switch

1. 1. Product Introduction

The curing light manufactured by Guilin Refi ne Medical Instrument Co., Ltd. uses the principle of light radiation to irradiate photosensitive resin in order to quickly cure it. This Curing Light is intended for use by trained dental professionals for the purpose of light curing dental resins and composites.

Patient target groups: adults and pediatrics.

Scope of application: suitable for dentistry, used to cure light-curing materials.

Product composition: main unit, LED light, light hood, charging base, battery and power adapter, etc.

The features of this product are:

1.1. A variety of different working modes and working hours can be set.

1.2. Constant light power output, the decrease of battery power will not affect the curing effect.

1.3. Large-capacity battery. While full charged, under 10 seconds of light per time, it can be used continuously for more than 500 times.


1.4. Curing Light Lenses:

Microcure Lens: creates a 2.5 mm diameter curing area to aid in position and tack curing.

Microcure Ball Lens: holds matrix against contact while curing composite for Class II restorations.

Endo Guide Lens: For apical retrofills, narrow preps.

1.5. Product schematic:

2. Replacement Instructions for Accessories

Accessories	Period of Replacement	Method of Replacement
Power adapter	Damaged	Replace directly according to the requirements of the manual.
Charging base	Damaged	Replace directly according to the requirements of the manual.
Battery	Damaged	Twist the cover at the bottom of the handpiece to remove and replace the battery. If the device starts up and operates normally after the battery is replaced, the replacement has been successfully completed.
Light hood	Damaged	Replace directly according to the requirements of the manual.

3. Technical paramaters

3.1. Technical Parameters of each model (table 1)

Dimensions	Φ86mm*223mm
Net weight	265g
Classified by power supply	Powered by rechargeable lithium battery
Charger	Input: 5V DC/1A Working: 3.7V lithium battery(ICR 18650)
Rechargeable lithium battery	Battery model: ICR 18650 Battery voltage and capacity: 3.7V/2600mAh The battery has overvoltage, overcurrent and short circuit protection.
Power supply (Adapter)	Input : 100-240V AC 50Hz/60Hz Output: DC 5V/1A Built-in fuse: T1A 250V
Performance of led lamp	a) 10W high power blue and purple LED lamp. b) Wave length: 385nm-515nm c) Category: class I.d)Shoot limit(AEL): 3.9x10 ⁻³ J
Light Intensity	1000mW/cm²-2500mW/cm²
Inspection method of LED lamp	When operating correctly, the LED lamp lights up, indicating that the lamp is in good condition. Dental resin materials commonly used in clinic can match the wavelength of this curing light, such as 3M, Dentsply and other resins.
Modes and time setting	Turbo mode: Blue and purple light, 2300-2500mW/cm2 for 1s, 3s Normal mode: Blue and purple light, 1000-1200mW/cm2 for 5s, 10s, 15s, 20s Ortho mode: Blue and purple light, 2300-2500mW/cm2 for 3s to work 10 times every cycle, or 5s to work 10 times every cycle, interval 1s Check mode: purple light for 30s, 60s
Effective area of irradiation light	75.4mm²

3.2. Working Environment

3.2.1. Operation: +5°C - +40°C

3.2.2. Relative humidity: 30% - 75%

3.2.3. Atmospheric pressure: 70kPa-106kPa

3.3. Equipment safety classification

3.3.1. Classified by type of protection against electric shock: Class II

3.3.2. Classified by degree of protection against electric shock: Type B application part

- 3.3.3. Classified according to the degree of protection against incoming liquid: Ordinary equipment (IPX0), non-waterproof
- 3.3.4. Classified by operation mode: Non-continuous operation, the maximum working time 100 seconds, the minimum deactivation time 10 minutes.
- 3.3.5. Classification of safety degree when using flammable anesthetic gas mixed with air or oxygen or nitrous oxide: Equipment that cannot be used with flammable anesthetic gas mixed with air or oxygen or nitrous oxide.

4. Installation and Demounting

Put on the disposable curing light sleeve to cover the entire curing light for cross infection control.
Clamp and install the light hood at the front of the curing light.

5. Function and Operation

- 5.1. When the machine is under stand-by mode, press the On/off key to light up the screen.
- 5.2. Mode setting: press the time/mode key for more than 2 seconds to select the mode; and every time the key is pressed, there's a beep. Without operating within 2 minutes, the screen will turn off and the main unit enters the sleep state.
- 5.3. Time setting: press the time/mode key to select time, see table 1 for settable curing time.
- 5.4. Point the light tip at the right position, press the switch button and there's a beep. The curing light will start to work as the mode setting. The display screen starts counting down from the selected time, and after the timing ends, the display screen displays the setting time.
- 5.5. Before the end of setting time, at any time pressing the on/off key can finish the work.
- 5.6. After the end of an irradiation cycle, you can immediately press the switch to start the next irradiation cycle. If the handle starts to become hot, please stop working until the handle is completely cooled down. It is recommended that the number of consecutive irradiations should not exceed 10 times.
- 5.7. The effective light intensity of this machine is many times higher than that of the halogen lamp curing light, and the curing depth of the light-curing composite resin is not less than 4mm for 10 seconds.
- 5.8. When the device is used continuously for a long time, the maximum temperature of the applied part may reach 49.2°C. After 100 seconds of continuous use, a 10-minute interval must be applied, or stop using when the temperature rises significantly.

6. Scaling function

- 6.1. Steam sterilization of this equipment is not permitted. Doing so may result in device damage. The handpiece can be handled in the same gesture as a pen in hand.
- 6.2. Before using this equipment, users need to put disposable isolation sleeve which is supplied by our company on the light guide of Curing Light to avoid contact between the main unit or other components and the patient's skin or oral mucosa and avoid cross-infection. After use, remove the disposable isolation sleeve and dispose of it in accordance with applicable legal requirements in your country as well as the hygiene regulations of the hospital or clinic.
- 6.3. The disposable isolation sleeve purchased by the user should meet the requirements of medical device regulations.
- 6.4. The accessory of the equipment should be cleaned by clean water or neutral sterilized liquid. Do not soak. Do not use highly volatile and diffuent solvent to clean this equipment, which can cause the signs on the control panel to fade.
- 6.5. For cleaning and disinfection, use disinfectant ethanol or disinfectant isopropanol that does not contain any additives (any quaternary ammonium salts). Use of other disinfectants may cause discoloration or cracking.
 - For details, contact the manufacturer of the disinfectant.
 - In this operation manual, disinfectant ethanol or disinfectant isopropanol that does not contain any additives (any quaternary ammonium salts) is called by "disinfectant alcohol".
- 6.6. Disposable isolation sleeves are not allowed to be reused to prevent cross-infection.
- 6.7. Please clean the resin remained on the top of the main unit after using to avoid affecting the lifespan or solidified effect.

7. Transportation, storage and maintenance

- 7.1. Transportation
 - 7.1.1. During transportation, avoid shocking and vibration, handle with care and avoid inversion.
 - 7.1.2. Do not mix with dangerous goods during transportation.
 - 7.1.3. Avoid sunlight, rain or snow during transportation
- 7.2. Storage
 - 7.2.1. Do not mix with toxic, corrosive, flammable and explosive materials during storage.
 - 7.2.2. The product should be stored in an environment where the relative humidity does not exceed 10% - 93%, the atmospheric pressure is 70kPa - 106kPa, and the temperature is -20°C - +55°C.
- 7.3. Maintenance
 - 7.3.1. This equipment should be handled with care, far away from the source of the earthquake, and should be installed or stored in a cool, dry and ventilated place.
 - 7.3.2. When the device is not in use, turn off the power switch and unplug the power plug; when not in use for a long time, charge the main unit every 3 months.
 - 7.3.3. To clean or disinfect, wipe the surface with clean water or disinfectant, do not soak.
 - 7.3.4. After use, check if there is any resin remaining on the surface of lamp, so as not to affect the life of the front joint or affect the curing effect.

8. Troubleshooting

Faulty	Possible Cause	Solutions
No reaction	1. Out of battery 2. Faulty of battery 3. Battery protection	1. Charge/Change a new battery 2. Change a new battery 3. Activate it by putting it on the base
Insufficient light intensity	There is resin remaining on the surface of the transparent cover of light tip.	1. Remove the resin 2. Change the transparent cover
Not charging when connected with adapter	1. The adapter is not connected well 2. Faulty of adapter or incompatible 3. The charging point is impurity.	1. Reconnect. 2. Change the adapter 3. Clean the inner connector of charging base by the alcohol
Less working life after fully charged	Capacity of battery decreased	Replace with a new battery
The mode indicator twinkles when charging	Low voltage	Back to normal after 15-minute charging

If the fault still cannot be resolved, please contact your local dealer or our company

9. EMC-Declaration of conformity

- 9.1. Instructions for use
 - The ME EQUIPMENT or ME SYSTEM use in hospitals or dental clinics.

- ⚠ Warning: Do not operate the device near active high-frequency (HF) surgical equipment or within the RF-shielded room of a magnetic resonance imaging (MRI) system, as strong electromagnetic disturbances may interfere with device function.
- ⚠ Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- ⚠ Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- ⚠ Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Note: The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

List of all cables

No.	Name	Length	Shielded or not	Detachable or not	Note
1	Power adapter output cable	1.2m	No	Yes	/

Replaceable accessories

No.	Name	Model Specification	Connection method	Note
1	Power adapter	5V/1A	Plug	/
2	Battery	ICR 18650	Plug	/

Essential performance

The curing light has neither life sustaining functions nor diagnostic of life supporting functions.

The following functions are observed:

- Curing of dental resins and composites
- Continuous light on handpiece

When the product is invalid or degraded due to electromagnetic disturbance, the user should stop using it immediately to ensure that there is no error caused by the product's performance failure or degradation. In this case, the user should remove the disturbance source or adjust the direction or position of the product, so that the product can work normally.

9.2. Technical description

- 9.2.1. Portable and mobile RF communications equipment may affect the performance of equipment, use of equipment should be avoided near strong electromagnetic interference, and do not use close to mobile phone, microwave oven, etc.
- 9.2.2. Use of equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- 9.2.3. Except for the cables sold by manufacturers of as spare parts of internal components, the use of accessories and cables other than those specified or provided by the manufacturer may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- 9.2.4. Use of accessories, transducers and cables other than those specified or provided by the manufacturer together with equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- 9.2.5. Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

Table 1

Guidance and Manufacturer's Declaration – Electromagnetic Emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Applicable

Table 2

Guidance and Manufacturer's Declaration – Electromagnetic Immunity		
Immunity test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV power supply lines ±1 kV signal input/output 100 kHz repetition frequency	±2 kV power supply lines Not applicable 100 kHz repetition frequency
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode	±0.5 kV, ±1 kV differential mode Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle	0 % UT; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
NOTE UT refers to the a.c. mains voltage measured prior to the application of the test level.		

Table 3

Guidance and Declaration – Electromagnetic Immunity						
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test frequency (MHz)	Band (MHz)	Services	Modulation	IEC 60601-1-2 test level (V/m)	Compliance level(V/m)
	385	380 – 390	TETRA 400	Pulse modulation 18 Hz	27	27
	450	430 – 470	GMRS 460, FRS 460	FM± 5 kHz deviation 1 kHz sine	28	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	9	9
	745					
	780					
	810					
	870	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28	28
	930					
	1720	1700 – 1990	GSM 1800;CDMA 1900;GSM 1900; DECT;LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	28	28
	1845					
	1970					
	2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28	28
	5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9
	5500					
	5785					

Table 4

Guidance and Manufacturer's Declaration – Electromagnetic Immunity				
Radiated RF IEC61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	Test Frequency	Modulation	IEC 60601-1-2 Test Level (A/m)	Compliance level (A/m)
	30 kHz	CW	8	8
	134,2 kHz	Pulse modulation 2.1 kHz	65	65
	13,56 MHz	Pulse modulation 50 kHz	7,5	7,5

10. Environmental Protection

Components	Hazardous substance or element					
	(Pb)	(Hg)	(Cd)	(Cr6+)	(PBB)	(PBDE)
Main units	○	○	○	○	○	○
Charging base	○	○	○	○	○	○
Adaptor	○	○	○	○	○	○
Mechanical element (Including screws, nuts, joint rings, etc.)	○	○	○	○	○	○
○ : Indicates that the content of this toxic and hazardous substance in all homogeneous materials of this component is below the limit requirement in SJ/T-11363-2006 "Limit Requirements for Toxic and Hazardous Substances in Electronic Information Products".						
X : It means that the content of this toxic and hazardous substance in at least one of the homogeneous materials of the part exceeds the limit requirement of SJ/T11363-2006. (This product complies with EU RoHS environmental protection requirements; currently there is no mature technology in the world that can replace or reduce the lead content in electronic ceramics, optical glass, steel and copper alloys)						
According to the "Administrative Measures for Restricting the Use of Hazardous Substances in Electrical and Electronic Products", "Regulations for Recycling and Disposal of Waste Electrical and Electronic Products" and relevant standards, please comply with the safety and use precautions of the products, and recycle or dispose of the products in appropriate ways according to local laws and regulations after the use of the products.						

11. After sales service

- 11.1. If this equipment cannot work normally due to quality problems from the date of sale, please show the warranty card and our company will be responsible for the maintenance. Please refer to the product warranty card for the warranty period and warranty scope.
- 11.2. This product does not contain replaceable parts. All maintenance, adjustment, calibration and modification of technical parameters of the product can only be carried out by our company's technicians or special repair shops. If customers need to repair themselves, our company can provide them on request Circuit diagrams, component lists, legends, calibration specifications, or other materials needed to help qualified technical personnel repair the equipment parts designated by our company, but we will not bear any consequences.
- 11.3. The user must use only original parts. For purchases, please contact the local distributor or the manufacturer. Do not use accessories from other brands, as this may cause damage to the curing light or result in other hazards.
- 11.4. After the power adapter, charging base and other accessories are damaged, users should not repair by themselves. Please buy new parts and replace them by yourself. If you need relevant information, please contact the manufacturer.

Note: We reserve the right to modify the equipment design, technical specifications, fittings, instruction manual, and original packing list at any time without prior notice. In the event of discrepancies between the blueprint and the actual equipment, the actual equipment shall be considered the standard.

Refer to the product label for the date of manufacture.

 **Guilin Refine Medical Instrument Co., Ltd.**
No.8-3, Information Industrial Park, High-Tech Zone, Qixing District, 541004 Guilin, Guangxi, PEOPLE'S REPUBLIC OF CHINA
Tel: +86-773-7796686 Email:refine@refine-med.com
Website: http://www.refine-med.com

 MedNet EC-REP C llb GmbH
Borkstrasse 10, 48163 Münster, Germany

Scan and log in to the website for more information.

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MedNet EC-REP C Ilb GmbH