SDI LUNA / LUNA FLOW / FLOW LV

ENG UNIVERSAL FLOWABLE COMPOSITE
BUL УНИВЕРСАЛЕН ТЕЧЕН КОМПОЗИТ
CES UNIVERZÁLNÍ TEKUCÍ KOMPOZIT
DAN UNIVERSAL FLOW KOMPOSIT
DEU UNIVERSAL FLIESSFÄHIGES KOMPOSIT
EL ΓΕΝΙΚΗΣ ΧΡΗΣΗΣ ΡΕΥΣΤΟ ΣΥΝΘΕΤΗ ΡΗΤΙΝΗ
FIN UNIVERSAALI JUOKSEVA YHDISTELMÄMUOVI
FRA COMPOSITE FLUIDE UNIVERSEL
HUN UNIVERZÁLIS FOLYÉKONY KOMPOZIT
ITA COMPOSITO FLUIDO UNIVERSALE
LIT UNIVERSALUS TEKANTIS KOMPOZITAS
NOR UNIVERSAL FLYTENDE KOMPOSITT
POL UNIWERSALNY KOMPOZYT TYPU "FLOW"
POR RESINA COMPOSTA UNIVERSAL FLUIDA
RON MATERIAL COMPOZIT FLUID UNIVERSAL
SLK UNIVERZÁLNY TEKUTÝ KOMPOZIT
SPA COMPOSITE FLUIDO UNIVERSAL
SWE UNIVERSAL LÄGVISKÖS KOMPOSIT

Made in Australia by SDI Limited
3-15 Brunsdon Street, Bayswater, Victoria 3153, Australia

Australia 1 800 337 003
Austria 00800 0225 5734
Brazil 0800 770 1735
France 00800 0225 5734
Germany 0800 1005759
Italy 00800 0225 5734
Spain 00800 0225 5734
UK 00800 0225 5734
USA & Canada 1 800 228 5166

EU REP
SDI Germany GmbH
Hansestrasse 85
51149 Cologne, Germany
T: +49 2203 9255 0

Last Revise: 06-2025

PAK41384 E

-25°C -77°F / 2°C 35°F / Allergies & Precautions / CE 0297 / Rx only / MD / Paper

ENG INSTRUCTIONS FOR USE

INTENDED USE
Luna Flow /Luna Flow LV
is designed to restore teeth
by dental professionals.
Luna Flow / Luna Flow LV flowable composite, is a versatile, radiopaque, flowable, light cured material. The composite is directly extruded into the cavity preparation for maximising adaptation to the preparation.

COMPOSITION OF LUNA FLOW:
Dimethacrylates (43% wt) (68% vol)
• Urethane dimethacrylate
• Triethylene Glycol Dimethacrylate
• Ethoxylated Bisphenol A Dimethacrylate
Fillers (56%wt) (32% vol)
(0.02 – 1 micron)
• Barium Glass
• Silicon Dioxide (hydrophobic fumed silica)
• Ytterbium Trifluoride
Initiators (< 1%)
Stabilisers (< 1%)
Pigments (< 1%)

COMPOSITION OF LUNA FLOW LV:
Dimethacrylates (47% wt) (72% vol)
• Urethane dimethacrylate
• Triethylene Glycol Dimethacrylate
• Ethoxylated Bisphenol A Dimethacrylate
Fillers (52%wt) (28% vol)
(0.02 – 1 micron)
• Barium Glass
• Silicon Dioxide (hydrophobic fumed silica)
• Ytterbium Trifluoride
Initiators (< 1%)
Stabilisers (< 1%)
Pigments (< 1%)

Luna Flow / Luna Flow LV composite is radiopaque with about 265% and 250% aluminium respectively (100% aluminium has a radiopacity equivalent to that of dentin and 200% aluminium is equivalent to enamel).

INDICATIONS FOR USE:
- Conservative Class I (including occlusal surfaces), II, III, IV restorations
- Class V restorations (cervical abrasion, root erosion, abfraction)
- Repair of defects in enamel and porcelain restorations
- Sealants (pit & fissure, implant)
- Radiopaque cavity liner under direct restorations
- Blocking out of undercuts
- Minor core build ups
- Repairing or sealing of temporary crowns
- Cementation of porcelain, ceramic, composite veneers
- Splinting (mobile teeth, fibre bridges)
- Cover stains.

CONTRA-INDICATIONS:
- **Do not** use if a dry field cannot be established
- **Do not** use for pulp capping.
- **Do not** use in conjunction with any eugenol-containing materials.
- **Do not** use in anyone who has an allergy to acrylics or who has an allergic reaction to the material of this product.

NOTE: This product may cause skin irritation to some people. In such cases, discontinue use and seek medical attention.

SIDE EFFECTS:
- Areas around the pulp should be protected with a suitable pulp/dentin protector to prevent any potential irritation (selectively apply a calcium hydroxide-based preparation in areas near the pulp).
- The components of Luna Flow / Luna Flow LV could, on rare occasions, cause sensitivity. In these cases, Luna Flow / Luna Flow LV should not be used.

INTERACTIONS:
- Eugenol containing substances may inhibit the polymerisation process of composite materials. Therefore, using such materials in conjunction with Luna Flow / Luna Flow LV must be avoided.
- The use of chlorhexidine, cationic mouthwashes, and plaque-disclosing substances may result in discoloration.

INSTRUCTIONS FOR USE:
Preparation recommendations:
1. Clean and isolate tooth employing standard techniques and instruments. Rubber dam is the preferred mode of isolation.
2. A conservative cavity should be prepared, employing standard techniques and instruments, to form a slightly undercut internal form. Margins of the cavity preparation should end in sound and supported enamel with no bevels in stress bearing locations. If bevelling is desired in a low stress location it should extend no further than 1mm at an angle of no greater than 45°. Pre-wedging is also recommended to ensure that the restored tooth will have an adequate contact point.
3. Prophy all surfaces to be etched including surfaces adjacent to the cavity with an oil free non-fluoride containing paste or a slurry of pumice and water.
4. Rinse thoroughly with water.

Shade Selection
5. Ensure the tooth surface is clean and dry before colour matching. Do not desiccate.
6. Use a VITA® shade guide to select your shade before (35-77°F) it rubber dam.

NOTE: Do not use uncured paste for colour matching as there is a slight colour change upon curing.

Recommendation:
When the adjacent teeth are significantly different in shade, try to match the lighter shade. For visual confirmation of shade selected, place small amounts of composite on a surrounding tooth surface and light cure before bonding.

Isolation
Isolation techniques must be used to prevent contamination. Rubber dam is the preferred mode of isolation. Do not use in a dry field cannot be established.

Direct Restoration
7. Thoroughly dry the surface to be etched with dry, oil-free air. Etch tooth surface using the 'selective etch' or 'total etch' technique with 37% phosphoric acid according to the manufacturer's instructions, or proceed to the following steps if using the 'self etch' technique.
8. Wash thoroughly with water.
9. Remove excess water. Keep moist. Avoid contamination e.g. saliva.
10. Apply bonding agent, such as Zipbond (SDI Limited), according to the manufacturer's instructions to saturate all internal surfaces.
11. Directly extrude the flowable composite in increments of 2mm or less.
NOTE: ensure that the dispensing tip is correctly attached to the syringe by twisting the tip hub securely onto the syringe.
NOTE: For best results, apply restorative at room temperature (23 ± 1°C, 73 ± 2°F).
⚠ **CAUTION:** Air bubbles may be present when extruding the material. In these instances, use a probe to burst air bubbles before light curing.
12. Light cure the composite using a suitable high power LED curing light (440-480nm wavelength with a minimum light intensity of 1,000 mW/ cm²) such as Radii-Cal CX (SDI Limited) for a minimum of 20 seconds in increments of 2mm.
13. Finish the restoration using standard techniques. Instruct patient not to eat or drink anything for at least one hour after the procedure.
NOTE: Cooling with water spray and effective suction whenever possible when finishing the restoration is recommended.

Additional notes:
Luna Flow / Luna Flow LV syringe tips are for single use only. Reuse will result in cross contamination.

CEMENTATION OF PORCELAIN, CERAMIC AND COMPOSITE VENEERS
Prior to cementation with Luna Flow / Luna Flow LV, the restoration must be pre-treated according to the manufacturer's instructions. Luna Flow / Luna Flow LV is only suitable for the cementation of restorations that have translucency and are of appropriate thickness to allow polymerisation light to pass through, to ensure that the light reaches Luna Flow / Luna Flow LV, so that it can cure completely.

STORAGE AND HANDLING
- Avoid excess heat.
- Store at temperatures between 2° and 25°C (35° - 77°F).
- For optimum freshness, refrigeration is optional.
- Product temperature should be between 20° and 25°C (68° and 77°F) when used.

SIDE EFFECTS:
- Areas around the pulp should be protected with a suitable pulp/dentin protector to prevent any potential irritation (selectively apply a calcium hydroxide-based preparation in areas near the pulp).
- Patients with known material sensitivity or allergic reactions to this product should not use it.
- For professional use only.
- Keep out of reach of children and pets.
- Do not take internally.
- Do not use after expiry date.
- Do not expose material to direct light.

⚠ **CAUTION:** Federal Law restricts this device to sale by or on the order of a dentist.

FIRST AID
Eye (contact): Wash thoroughly with water and seek medical attention.
Skin (contact): Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.
Ingestion: Rinse mouth with water and discard. Drink plenty of water/ milk. Do NOT induce vomiting. Seek medical attention if symptoms persist.
Inhalation: No symptoms expected under normal conditions of use. If affected, move to well ventilated area.

If a serious incident has occurred due to this product, please report this incident to the manufacturer of this product and to the competent authority in your state or country. Dispose in accordance with local regulations.

SDS available at www.sdi.com.au or contact your regional representative.

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