



pro•rinse™




- Improved debris removal- cannula with side port opening
- Sterilized probe and rounded-end for safe use
- With Luer-Lock system for convenient use

**+
WE
KNOW
ENDO.**

DENTSPLY
MAILLEFER


PRORINSE™ is an irrigating probe with a side port dispersal and a closed, rounded-end for root canal and periodontal pocket.



Periodontal Pockets

Side port produces flushing motion to irrigate pockets

Rounded tip protects tissue










Root Canals

Design produces upward flushing motion for canal irrigation

Side port dispersal reduces risks of debris from being expressed through the apex compared to open-ended probes

Closed, rounded end reduces risk of apex damage compared to open-ended probes



	Box of 40 Probes	Box of 100 Probes	100 Probes + 100 Syringes	Outer Diameter	Gauge
	A082100023400	A082100023100	A082100023S00	0.63 mm	23
	A082100024400	A082100024100	A082100024S00	0.56 mm	24
	A082100025400	A082100025100	A082100025S00	0.51 mm	25
	A082100028400	A082100028100	A082100028S00	0.36 mm	28
	A082100030400	A082100030100	A082100030S00	0.30 mm	30

Dentsply Nordic Baltic
 Box 2024, SE-128 21 Skarpnäck
 Tel +46 (0)8 6856505
 Fax +46 (0)8 4476055
 www.dentsply.eu

For local contacts, please visit
 www.dentsply.se
 www.dentsply.fi
 www.dentsply.no
 www.dentsply.eu
 www.dentsply.dk
 LV + 371 29 455 945
 LT +370 605 80858

